

## Spitfire EX2 3-Wheel Travel Scooter with 12AH Battery



INSET A



BLUE



SILVER




INSET D



INSET E



INSET F



INSET G



INSET H

- 6 sets of interchangeable color panels in red, blue, silver, black, white and pink
- Swivel mount cup holder included
- Extra bright LED headlight
- Charge ports located on the tiller and battery pack for easy charging
- Height adjustable swivel seat
- Adjustable padded armrests with reflectors
- Large plastic basket with removable liner
- Disassembles into 5 easy to carry pieces
- Delta tiller for easy throttle control

PRODUCTS	Description
SFEX2317FS-12	KO800, 1/ea
SFEX2317FS-21	KO800, 1/ea

SPECIFICATIONS	Description
Weight Capacity	300 lbs.
Max Speed	5 mph
Maximum Range	9 Miles   15 Miles
Max. Climbing Angle	6°
Turning Radius	46"
Ground Clearance	3"
Overall Dimensions	42" x 20" x 36"
Seat Width	17"
Seat Depth	16"
Back Height	16"
Seat-to-Floor Height	19" – 21"
Motor	24V x 270W x 4700 rpm
Controller	45A
Batteries	12AH   21AH
Battery Charger	2A Offboard
Brakes	Electromagnetic
Freewheel Mode	Yes
Front Wheels	7.5" x 2"
Rear Wheels	8.5" x 3"
Base Weight	60 lbs.
Battery Weight	22 lbs. (12AH)   28 lbs. (21AH)
Seat Weight	17 lbs.
Total Weight	99 lbs. (12AH)   105 lbs. (21AH)

**SPECIFICATIONS**

**Description**

Warranty on Frame	Limited Lifetime
Warranty on Electronics	24 Months
Warranty on Batteries	12 Months

**ACCESSORIES**

**Description**

SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea