

XLR 14



The XLR 14 Power Chair is equipped with large 14" Drive Wheels. The self-adjusting Mid-Wheel Drive suspension system is matched with four independently suspended front and rear casters. This provides the added stability of having six points of contact with the ground. When accompanied with the full independent suspension system, the XLR14 provides a smooth driving experience that will handle virtually any driving surface. A natural fit with our power tilt seating system, which offers a 50° of tilt. The XLR 14 also has a flip forward seat for easy battery access to our Xtra Long Range battery options, and comfortably accommodates a 300 lb. weight capacity.

- + Xtra Long Range Battery Capacity
- + 300lb Weight Capacity
- + 14" Mid-wheel Drive w/ Suspension System
- + Independently Suspended Front & Rear Casters
- + Power Tilt Seating System
- + Manual Semi Recline Back



- Red
- Blue
- Pewter
- Black

Model Name: XLR14-CS
Medicare Billing Code⁴: K0836

SPECIFICATIONS

WEIGHT CAPACITY	300 lbs.
WIDTH¹	24.8"
LENGTH¹	43"
UNIT WEIGHT³	283 lbs.
MAX SPEED	4 mph

GROUND CLEARANCE	2.75"
SAFE CLIMBING ANGLE	6° 10%
FRONT REAR WHEELS	6" Solid 5" Solid
DRIVE WHEEL	14" Solid
SUSPENSION TYPE	Active/Independent

BATTERY OPTIONS	Group 34	Group 24
BATTERY WEIGHTS	39 lbs.	51 lbs.
RANGE² UP TO	16 miles	20 miles
CHARGER	8A Off-board	
STFH	21" - 25"	

¹Due to manufacturing tolerances and continued improvement, this specification can be subject to variance. Dimensions are for the power base only, overall dimensions can change depending on seating and accessory selection. ²Varies with user weight, terrain, battery type, battery condition and tire pressure. ³Unit weight includes batteries and seating system. ⁴HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.