





6Runner 10

The 6Runner 10 is a mid-size power chair that is design for the daily user. It sports a mid-wheel drive system and has the added stability of having six points of contact with the ground. The 6Runner 10 also has our stable four post seating system and suspension system for smooth drive operation. The unit has easy battery access, an 18in. wide reclining seat and 300 lb. weight capacity.

- + 300lb Weight Capacity
- + Mid-wheel Drive
- + Full Suspension
- + Easy Battery Access
- + Four Post Seating System



Red

Model Number: 888WNLM Medicare Billing Code: K0823

HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.

Technical Specifications

0 110' ' (1)// 11)		40" 00" 40"
Overall Dimensions (L x W x H)	in	42" x 23" x 46"
Number, Size of Tires	Front, Rear	2, 7" (2, 10") 2, 7'
Suggested User Weight ¹	lb	Max. 300
Weight of Heaviest Part ²	lb	-
Weight ³	lb	191
Each Battery Module Weight	lb	23
Battery Size	Type x PCS	U1 x 2
Charger	Amp	Off board (5A)
Output Power of Motor	HP	.60 x 2
Maximum Speed	MPH	5
Maximum Climbing Angle/	%(°)	10(6)
Safe Climbing Angle ⁴		
Range (Per Charge) ⁶	Mile	10 - 11
Turning Radius	in	15.75"
Model #		888WNLM
Standard Seat	$W \times D \times BH$	SV191821B

Note: ¹ On level Road. ² Without basket and seating. ³ With basket, seat & batteries. ⁴ Base on user weight of 200LB. ⁵ Based on user weight of 165LB. ⁶ After the batteries and mechanical moving parts fully break in.