



Trident
Standard Power Wheelchair Front-Wheel Drive
 HCPCS: K0823

- Shroudless design offers easy access to internal components and batteries
- Semi-reclining seat with adjustable headrest offers three height adjustments and two depth options
- Armrests are padded, adjustable and removable
- Available with 18" or 20" width Captain's Seat
- Quick-release pins allow easy seat removal
- Adjustable, flip-up footplate includes rubber front bumper to prevent damage
- Positioning belt included
- Seat position can be adjusted forward and backward
- Programmable controller
- Optional swing-away controller arm sold separately
- Optional elevating leg rests sold separately
- Black, flat-free non-marking wheels
- Large, front anti-tip wheels feature suspension for safety, smoother ride and enhanced obstacle handling
- Controller is protected by metal frame
- Easy free-wheel operation
- Slide-out battery tray

PRODUCTS	Description
2850-18	18" Captain's Seat, 1/cs
2850-20	20" Captain's Seat, 1/cs

SPECIFICATIONS	Description
HCPCS	K0823
Weight Capacity	300 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	24"
Climbing Angle	7.5°
Ground Clearance	2"
Dimensions	42.25" (L) x 25.25" (W)
Floor-to-Seat Height	21" - 24"
Controller	PG VR2 50A
Motors (2)	24V x 320W x 4600 rpm
Batteries (pair)	12V x 36AH (U1)
Charger	4A Offboard
Brakes	Electromagnetic
Caster Wheels	8" x 2" Solid
Drive Wheels	14" x 3" Flat-free
Base Weight	82 lbs.

SPECIFICATIONS	Description
Battery Weight (pair)	53 lbs.
Seat Weight	51 lbs.
Total Weight	186 lbs.
Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)
Warranty of Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	14 Months
Warranty on Batteries	12 Months

ACCESSORIES	Description
SF8030	Cane/Crutch Holder, 1/ea
SF8010	Oxygen Tank Holder, 1/ea
TRID-2A37-LINX	Swing-Away Controller Arm, 1/ea
TRID-2A47	ELR Seat Bracket, 1/ea
SFELR	ELR Hangers, 1/ea
LK3JELR	ELR Set, 1/pr

PARTS	Description
TRID-2A01	1/ea
TRID-2A02	1/ea
TRID-2A39	1/ea
TRID-2A05	1/ea
TRID-2A43	1/ea
TRID-2A44	1/ea
TRID-2A38	1/ea
TRID-2A36	1/ea
S350169701	1/ea
TRID-2A20	1/ea
TRID-2A22	1/ea
TRID-2A23	1/ea
TRID-2A25	1/ea
TRID-2A24	1/ea
TRID-2A26	1/ea
TRID-2A32	1/ea
TRID-2A33	1/ea
LRA402208	1/ea
LRA402307	1/ea
TRID-2A45	1/ea
TRID-2A16	1/ea
TRID-2A17	1/ea
TRID-2A18	1/ea
TRID-2A19	1/ea
TRID-2A34	1/ea
TRID-2A13	1/ea
STDS1011	1/ea
TRID-2A14	1/ea
TRID-2A1	
TRID-2A10	1/ea
TRID-2A11	1/pr
TRID-2A12	1/ea
TRID-2A28	3/pkg
TRID-2A29	1/ea
TRID-2A30	1/ea
TRID-2A31	1/pr
TRID-2A41	1/ea
TRID-2A42	1/ea