



Cirrus Plus EC
Folding Power Wheelchair Rear-Wheel Drive

HCPCS: K0822



FIGURE A

- Foldable, lightweight frame for easy portability (Figure A)
- Attractive steel cross-brace frame
- Armrests are padded, adjustable and removable
- Adjustable-length controller mount
- Adjustable-tension back upholstery
- Swing-away footrests with heel loops
- Optional elevating leg rests (Item #CPELR) sold separately
- Positioning belt included
- Padded back upholstery with storage pouch and cushioned seat
- Programmable controller
- Flat-free, non-marking tires
- Easy free-wheel operation
- Manual wheel locks
- General foam seat cushion included
- Anti-tip wheels

PRODUCTS	Description
CPN16FBA	CirrusPlusECPowerChair 16FlpBkArm, 1/cs
CPN18FBA	CirrusPlusPowerChair 18in FlpBkArm, 1/cs
CPN20FBA	CirrusPlusPowerChair 20FlpBkArm, 1/cs

SPECIFICATIONS	Description
HCPCS	K0822
Weight Capacity	300 lbs.
Top Speed	5 mph
Maximum Range	15 miles
Turning Radius (16" 18" 20")	30" 32" 33"
Climbing Angle	6°
Ground Clearance	4"
CPN16FBA Dimensions	42.5" (L) x 24" (W)
CPN18FBA Dimensions	42.5" (L) x 26" (W)
CPN20FBA Dimensions	42.5" (L) x 28" (W)
Floor-to-Seat Height	18.5"
CPN16FBA Seat Dimensions	16" (W) x 16" (D)
CPN18FBA Seat Dimensions	18" (W) x 16" (D)
CPN20FBA Seat Dimensions	20" (W) x 16" (D)
Controller	PG VSi 50A
Motors (2)	24V x 320W x 4600 rpm
Batteries (pair)	12V x 36AH (U1)

SPECIFICATIONS

Description

Charger	4A Offboard
Brakes	Electromagnetic
Caster Wheels	8" x 2" Solid
Drive Wheels	12.5" x 2.25" Flat-free
Base Weight	90 lbs.
Battery Weight (pair)	56 lbs.
Total Weight	146 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller	24 Months
Warranty on Drive Train Components	24 Months
Warranty on Batteries	12 Months