



Titan
Transportable Powerchair

HCPCS: K0821









INSET A

INSET B

INSET C

INSET D

- Height-adjustable swivel seat features fold-down backrest (Figure A)
- Includes interchangeable color panels in red and blue
- Seat back reclines and features adjustable headrest
- Adjustable-length controller mount (Figure B)
- Armrests are padded, adjustable and removable
- Available with 16", 18", 20" and 22" wide Captain's Seat
- Optional swing-away controller arm (Item #AA4100) sold separately
- Solid, flip-up footplate with adjustable height, depth and angle
- Positioning belt included
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately
- Programmable controller
- Flat-free, non-marking tires (Figure C)
- Easy free wheel operation (Figure D)
- Anti-tip wheels

·	
PRODUCTS	Description
TITANLB18CS	Titan K0821 FWD 18 x 18 Capt.Seat, 1/cs
TITAN18CS	Titan K0821 FWD 18 x 18 Capt.Seat, 1/cs
TITAN20CS	Titan K0821 FWD 20 x 18 Capt.Seat, 1/cs
TITAN1616	Titan K0821 FWD 16 x 16 Capt.Seat, 1/cs
TITAN1618	Titan K0821 FWD 16 x 18 Capt.Seat, 1/cs
TITAN1816	Titan K0821 FWD 18 x 16 Capt.Seat, 1/cs
TITAN2020	Titan K0821 FWD 20 x 20 Capt.Seat, 1/cs
TITAN22CS	Titan K0821 FWD 22 x 20 Capt.Seat, 1/cs
SPECIFICATIONS	Description
HCPCS	K0821
Weight Capacity	300 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	25"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions	40" (L) × 24" (W)
Floor -to-Seat Height	20.5" - 23"
Controller	Dynamic Shark 40A
Motors (2)	24V x 200W x 3900 rpm
Batteries (pair)	12V x 35AH (U1)
Charger	3A Offboard



SPECIFICATIONS	Description
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Drive Wheels	10" x 3" Flat-free
Base Weight	68 lbs.
Battery Weight (pair)	44 lbs.
Seat Weight	43 lbs.
Heaviest Piece Disassembled	55 lbs.
Total Weight	158 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronic Controller/Driv	e Train Components 14 months
Warranty on Batteries	12 Months
ACCESSORIES	Description
AA4100	Swing-Away Joystick, 1/ea
AE2500	ELR Bracket, 1/ea
LK3JELR	ELR Set, 1/pr
AB1000	Arm Bag, 1/ea
AZ0060	Cup Holder, 1/ea