

Titan LTE Portable Powerchair



INSET A



INSET B



INSET C




INSET D



INSET E



INSET F



INSET G

- Disassembles for transport & storage
- Flip-back armrests with height, width and angle adjustments
- Interchangeable color panels in red & blue
- 300 lbs. weight capacity
- Flat-free, non-marking tires
- Bluetooth electronics for easy troubleshooting and diagnosing

PRODUCTS

Description

TITANLTE-18FS	Titan LTE Portable Powerchair, K0821, 1/ea
---------------	--

SPECIFICATIONS

Description

Product Code/SKU	TITANLTE-18FS
Product Type	Power Wheelchair
Seat Width	18"
Seat Depth	18"
Seat to Floor Height	20" - 22"
Turning Radius	27"
Steering Type	Directional Joystick
Max Speed	4.2 mph
Anti Tip Wheels	Standard
Batteries	(2) 12V x 20AH
Freewheel Mode	Yes
Brakes	Electromagnetic
Casters	8" x 2"
Charger	3A Offboard
Climbing Angle	6°
Controller	Dynamic Linx 40A
Maximum Range	8 miles
Motor	200W x 3900 rpm
Drive Wheels	9" x 3"
Ground Clearance	2.5"
Heaviest Piece	39 lbs.
Base Weight (w/out battery)	70 lbs.
Seat Weight	30 lbs.
Battery Weight	22 lbs.
Base Weight	92 lbs.
Warranty	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries

ACCESSORIES

	Description
SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea

PARTS

	Description
LINX-JOYSTICK	Joystick LiNX Electronics, 1/ea
DC GPL55470	Joystick Tip, 1/ea
GLM-BUS100-A	Linx Joystick Cable, 1/ea
LDM400001-07	Joy Stick Bracket Knob, 1/ea
LDM400001-04	Joystick Arm, 1/ea
LDM400001-01	Joystick Screw, 1/pr
LDM400001-05	Joystick Bracket sub, 1/ea
DLX-PM40-A-LTE	Controller, 1/ea
GLM-EXT064-A	Cntrlr To Shrd Cbl, 1/ea
C 10-085-00500	Controller Cable Clamp, 1/ea
DST-013	SEAT W/SEAT PLATE W/O ARMRESTS, 1/ea
LDR410011L	Armrest Assembly Left, 1/ea
LDR410011R	Armrest Assembly Right, 1/ea
LDR410110	Arm Pad, 1/ea
C07-069-00201	SEAT POST, 1/ea
TITANLTE-11	Seat Post Bolt, 1/ea
DST-016	SEAT PLATE, 1/ea
C09-085-00200	Motor, Gearbox, Brake Left, 1/ea
C09-085-00100	Motor, Gearbox, Brake Right, 1/ea
C09-085-00700	Brake Only Left, 1/ea
C09-085-00800	Brake Only Right, 1/ea
C88-085-00100	Motor Key Way, 1/ea
531132B-HD	Rear Wheel, 1/ea
TITANLTE-04	Rear Wheel Mntg Bolt & Washers, 1/ea
TITANLTE-01	Articulating Bar, 1/ea
TITANLTE-13	Articulating Bar Mounting Hdwr, 1/ea
C26-084-00100	Dust Cover, 1/ea
DST-010	FORK STEM BEARING SET, 1/ea
CO3-017-00800	Fork, 1/ea
CO5-085-00100	Caster, 1/ea
SEC9031-B	Caster Bearing, 1/ea
TITANLTE-03	Caster Axle, 1/ea
C08-085-00100	Front Shroud, 1/ea
C08-085-00500	Rear Shroud, 1/ea
TITANLTE-16	Handle, 1/ea
TITANLTE-07	Red Decorations Set, 1/ea
TITANLTE-08	Blue Decoration Set, 1/ea
VENT-86	Battery Tray/Shroud Knob, 1/ea
C01-085-00100	Front Frame, 1/ea
C01-085-00200	Rear Frame, 1/ea
TITANLTE-09	Footrest & Mount, 1/ea
TITANLTE-14	Contact Block Assy & Hardware, 1/ea
TITANLTE-06	Anti-tip Assembly, 1/ea
TITANLTE-05	Frame Lock Bolt Assembly, 1/ea
LRM412123	20 AH BATTERY (TWO NEEDED), 1/ea
C 17-085-00100	Battery Bag, 1/ea
C31-005-00400	Battery Strap, 1/ea
TITANLTE-12	Wire Harness, 1/ea
LRM302107	20 Amp Circuit Breaker, 1/ea