

Vision CF

P322A

Front-Wheel Drive

Powerbase Wheelchair



Vision CF

P322A

Front-Wheel Drive

Powerbase Wheelchair

- ◆ Weight capacity 300 lbs
- ◆ Detachable, height and width adjustable armrests
- ◆ Removable, replaceable seat and back covers
- ◆ 18" seat width, ergonomic seating design
- ◆ Semi-reclining back rest 90°-135°
- ◆ Simplified frame construction
- ◆ Compact, lightweight design
- ◆ Easy access to components
- ◆ 9" Drive wheels, 6" castors
- ◆ Excellent maneuverability
- ◆ Flip up footplate



Technical Specs

| | |
|------------------|---|
| Max Speed | up to 4 mph |
| Range | up to 18 miles |
| Turning Radius | 21" |
| Ground Clearance | 2.4" |
| Motor | DC 24V, 130W |
| Controller | Dynamic LiNX |
| Battery | 12V/U1(34AH) x 2pcs |
| Charger | 4A off-board |
| Gradient | 6° |
| Casters | 6" foam filled tires |
| Drive Wheels | 9" foam filled tires |
| Brakes | Intelligent, regenerative, electro-magnetic |

Dimensions

| | |
|-------------------------|-----------|
| Seat Width | 18" |
| Seat Depth | 16" |
| Seat Height | 20" - 22" |
| Backrest Height | 17" |
| Min. Footrest Extension | 16.5" |
| Max. Footrest Extension | 18.5" |
| Overall Length | 33" |
| Overall Height | 42" |
| Overall Width | 24" |
| Total Weight | 154 lbs |
| Capacity | 300 lbs |

www.meritsusa.com



730 N.E. 19th Place, Cape Coral, FL 33909

TEL: 1-900-963-7487, 1-239-772-0579

Sales fax: 1-239-574-2661