



Seat Colors  
Charcoal/Gray Vinyl

Body Colors  
Red  
Blue

**HD Version  
450 lb.  
weight capacity!**

# Roll

golden  
**COMPASS™  
HD**

The heavy-duty version of the original Compass offers users a robust 450 lb. weight capacity with all of the same great features including the patented springless articulating chassis and true center-wheel drive design for unmatched maneuverability in tight spaces. The Compass HD comes in a variety of seat sizes up to 24" wide.




**YOUR LIFE IN MOTION**

**GOLDEN**  
technologies

## Golden Compass™ HD Specifications

Model	GP620CC	GP620SS
Model Name	Compass HD	Compass HD
Medicare Code	K0825	K0824
Weight Capacity	450 lbs	450 lbs
Drive Wheels	Center	Center
Maximum Speed	4 mph	4 mph
Operating Range <sup>1</sup>	12 Miles	12 Miles
Ground Clearance	3.5"	3.5"
Turning Radius	20"	20"
Type Batteries	2-NF22	2-NF22
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote
<b>Dimensions:</b>		
Length	40"	40"
Width	25"	25"
Height (ground to top of back of seat)	38.25" - 42.75"	38.25" - 42.75"
Height (ground to top of back of headrest)	47.25" - 53.75"	47.25" - 53.75"
Ground to Top of Seat	21.25" - 24.25"	19.75" - 22.25"
Number of Seat Height Positions	5	5
Size of Increments	0.5"	0.5"
<b>Weight of Unit: (assembled)</b>		
Base	101 lbs	101 lbs
Seat w/Arms	51 lbs	48 lbs
Batteries	37 lbs ea	37 lbs ea
<b>Tires:</b>		
Center	10"	10"
Casters	6"	6"
<b>Standard Seat:</b>		
	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 34.5"
Width x Depth	20" X 20", 22" X 20", 24" X 20"	20" X 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
<b>Optional Seats:</b>		
	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 34.5"
Width	18", 20", 22" or 24"	18", 20", 22" or 24"
Depth	18" or 20"	16" or 18"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
<b>Optional Seats 2:</b>		
	N/A	Pro-Flex Seat
Width	N/A	18", 20" or 22"
Depth	N/A	18" or 20"
Color	N/A	Black
<b>Footrest:</b>		
	Yes	Yes
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 - 35 degrees	0 - 35 degrees
<b>Warranty:</b>		
Frame	5 Years	5 Years
Drive Train	13 Months	13 Months
Electronics <sup>2</sup>	13 Months	13 Months

<sup>1</sup> Battery range will vary due to rider weight, drive surface and terrain.  
<sup>2</sup> Electronics warranty excludes batteries.



## Power Chair Accessories and Options

- Captain's seat (GP620CC)
- Pro Flex seat (GP620SS)
- Van/pan seat (GP620SS)
- Swing-away legrests
- Elevating legrests
- Stump support
- Footplate extender
- Cane/crutch holder kit
- Walker holder kit
- Oxygen tank holder kit
- Pack N Go

## 2 exciting colors



Red



Blue

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.