



USB  
Charger  
Port

LC-525i shown in  
Aria® Cool Grey fabric

MADE IN  
**USA**  
OF FOREIGN & DOMESTIC COMPONENTS

# INFINITY COLLECTION

The Infinity Collection from Pride® Power Lift Recliners features true infinite-position seating and pocketed coils designed to deliver comfort, style and performance. Available in four different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.

AMERICA'S #1  
POWER LIFT RECLINER

**POWER**  
LIFT  
RECLINERS

Protect your  
investment with  
the Pride® Lift  
Chair Care Kit



US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —

# INFINITY COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER LIFT  
RECLINERS



LC-525i shown in  
Ultraleather™ Pecan fabric



|                            | Petite Wide   | Small   | Medium  | Large   |
|----------------------------|---|---|---|---|
| <b>SPECIFICATIONS</b>      | <b>LC-525iPW<sup>2</sup></b>  | <b>LC-525iS<sup>2</sup></b>   | <b>LC-525iM<sup>2</sup></b>   | <b>LC-525iL<sup>2</sup></b>   |
|                            | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'2" & Below | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'3" & Below | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'4"- 6'0" | True Infinite-Position<br>Zero Gravity<br>Trendelenburg<br>Chaise Lounger<br>5'9"- 6'2" |
| <b>Weight Capacity</b>     | 375 lbs.  | 375 lbs.  | 375 lbs.  | 375 lbs.  |
| <b>Seat-to-Floor</b>       | 19"   | 19"   | 20"   | 21"   |
| <b>Seat Depth</b>          | 19"   | 19"   | 20"   | 22"   |
| <b>Seat Width</b>          | 24"   | 19"   | 20"   | 22"   |
| <b>Top of Back to Seat</b> | 28"   | 26"   | 29"   | 30"   |
| <b>Back Style</b>          | Blown Fiber<br>Split Back   | Blown Fiber<br>Split Back   | Blown Fiber<br>Split Back   | Blown Fiber<br>Split Back   |
| <b>Fully Padded Chaise</b> | Yes   | Yes   | Yes   | Yes   |
| <b>Heat and Massage</b>    | Optional  | Optional  | Optional  | Optional  |
| <b>Seating Type</b>        | Pocketed Coil<br>Springs  | Pocketed Coil<br>Springs  | Pocketed Coil<br>Springs  | Pocketed Coil<br>Springs  |
| <b>Distance From Wall</b>  | 19"   | 20"   | 20"   | 21"   |
| <b>Overall Chair Width</b> | 38"   | 34"   | 34.5"   | 37"   |
| <b>Chair Weight</b>        | 130 lbs.  | 130 lbs.  | 130 lbs.  | 130 lbs.  |
| <b>Warranty</b>            | 7-Year Prorated <sup>*</sup>  | 7-Year Prorated <sup>*</sup>  | 7-Year Prorated <sup>*</sup>  | 7-Year Prorated <sup>*</sup>  |

## INFINITY FABRIC SELECTION

QS

### DURASOFT



EMBER



OAT



DEEP SKY



TIMBER

### LEXIS STA-KLEEN® sta\*kleen. PERFORMANCE FABRIC



BLACK



MUSHROOM



BURGUNDY



CHESTNUT

### ULTRALEATHER™

UF  
ULTRALEATHER



BUFF



FUDGE



GARNET



PECAN



CHARCOAL

### ARIA<sup>3</sup> CRYPTON. super fabric



ESPRESSO



COOL GREY



RED



LAZULI



SAND

QS

Look for the Quick Ship  
Indicator for the shortest  
lead times available.

## BENEFITS

- Premium furniture grade laminate/hardwood frames for added strength and durability
- Pocketed coils deliver individualized support and relaxation
- Easy-to-use hand control with USB charger port *(not available when adding optional heat and massage)*
- Standard lumbar pillow
- Standard head and arm covers
- Large, dual pockets
- Easy to remove back
- Integrated battery backup in case of power failure<sup>1</sup>
- Transformer features self-diagnostic electronics with 9' cord

## OPTIONS\*

- Deluxe heat and massage
- Residential cleaning kit
- Head pillow
- Protective pad
- Lift chair overlay

**Additional lead time may apply for some options.**

*\*Additional charges apply for all options.*

1. Batteries not included.
2. Available with deluxe heat and massage only. Standard heat and massage not available.
3. Crypton<sup>®</sup> fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

<sup>\*</sup>Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.