


# Outlander Full

## Exterior Hitch Mount Platform Lift



Toggle Switch  
for ease of up  
and down use

Integrated License  
Plate Display

LED Puddle  
Light for ease of  
use at night

### DESIGNED FOR SMALL SCOOTERS AND POWER CHAIRS

The simple hitch-based installation requires little to no modification to your vehicle and allows for easy removal. The easy load/unload of an unoccupied mobility device onto the platform enables access from either end of the lift.

Consumers have the choice of a retractable strap kit, hold-down arm or hold-down bar.

# 25"

Max Travel Distance from Ground

# 350

Pound Capacity

# LOW

Maintenance

# 28"

Maximum Product Width

#### OUTLANDER FULL PRODUCT SPECIFICATIONS

Weight Capacity	350 lb (159 kg)
Vehicle Platform(s)	Full-Size Van, Minivan, Truck, Sedan, SUV and Crossover
Accommodates	Scooter, Power Chair
Securement Option(s)	Automatic Hold-Down Arm, Automatic Hold-Down Bar, Strap Kit
Platform Size	31.1" x 42" (79 cm x 106.7 cm) with a max. usable space of 30.3" x 36" (77 cm x 91.4 cm)
Power	12V DC
Hitch Class	III, IV
Max Travel Distance from Ground	25" (63.5 cm)
Option(s)	Battery Pack with Charger, Swing-Away Adaptor
Warranty	3-Year Transferable

harmar.com/paragon | 800.833.0478



7FEB20 | MKT-000143 REV A