

Travelite Transport Chair



FIGURE A



FIGURE B



FIGURE C



RETAIL BOX

- The Travelite Chair in a Bag is a lightweight yet sturdy aluminum transport chair in gunmetal gray
- It is easily transportable and comes with a compact carry bag with wheels making it ideal for traveling (Figure A - B)
- Folds for storage and transportation (Figure C)
- It is ideal for use indoors and outdoors
- 18" square nylon seat comes with lift-up straps
- Comes with wheel locks 8" solid tires
- Comes with padded, flip-up desk arms and height adjustable footrests

PRODUCTS

Description

TC005GY	Chair in a Bag, 1/cs
---------	----------------------

SPECIFICATIONS

Description

Overall Dimensions	22.5" (W) (open) 11.5" (W) (closed) x 37.5" (H)
Seat-To-Floor Height	19"
Back Height	17.5"
Seat Dimensions	18" (W) x 18" (D)
Overall Length (with riggings)	33"
Weight	20.5 lbs.
Carton Shipping Weight	23.2 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime