



Super Light, Folding Transport Chair with Carry Bag and Flip-Back Arms



FIGURE A

- Padded flip-back armrests
- Fold-down footrests
- Weighs only 19 lbs
- Easy to push or transport
- Aluminum frame with lightweight nylon upholstery
- Composite wheels are lightweight and maintenance free
- Folds like a conventional folding chair (Figure A)
- Comes with carry bag (item #835) with shoulder strap (Figure B)
- Standard rear wheel locks and seat belt

PRODUCTS

Description

SL18	Silver Frame and Blue Upholstery, SuperLiteTransportChair 181/cs, 1/cs
------	--

SPECIFICATIONS

Description

Overall Dimensions	23" (W) (open) x 37" (H)
Seat-To-Floor Height	20"
Back Height	17"
Seat Dimensions	18" (W) x 17" (D)
Overall Length	30"
Weight	19 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime

PARTS

Description

SL182A00	Front Caster, 1/ea
SL182A01	Rear Caster, 1/ea
SL182A02	Front Axle Bolt, 1/ea
SL182A03	Rear Axle Bolt, 1/ea
SL182A04	Footrest, 1/ea
SL182A05	Frame Dustcaps, 1/ea
SL182A06	Frame Spacer, 1/ea
SL182A07	Arm Pad, 1/ea
SL182A08	Seat Upholstery, 1/ea
SL182A09	Right Arm Mounting Bracket, 1/ea
SL182A0L	Left Brake, 1/ea
SL182A0R	Right Brake, 1/ea
SL182A10	Seat Belt, 1/ea
SL182A11	Folding Hinge, 1/ea

PARTS	Description
SL182A12	Support Bar, 1/ea
SL182A13	Caster Fork, 1/ea
SL182A14	Fork Bearing, 1/ea
SL182A15	Support Chain, 1/ea
SL182A16	Chain Clips, 1/ea
SL182A17	Caster Lock Knob, 1/ea
SL182A18	Back Upholstery, 1/ea
SL182A19	Left Arm Mounting Bracket, 1/ea
SL182A20	Footrest Hardware, 1/ea
SL182A21	Support Bar Hardware, 1/ea
SL182A22	Plastic Bolt Cover, 1/ea
SL182A23	Upholstery Screws and Washers, 1/ea
SL182A24	Hand Bar Pad, 1/ea
SL192A25	Rail Guide