

Fly-Lite Aluminum Transport Chair



FIGURE D

- Weighs only 16.2 lbs without footrests; 18.9 lbs with footrests
- Aluminum frame is lightweight and strong
- Comes with carry pocket on backrest
- Comes with seat belt for added safety
- Attractive plaid style nylon upholstery (Figure A)
- Comes with handy carry pocket on back (Figure B)
- Tool-free, height-adjustable swing-away footrests come standard
- Optional leg rests available (Model# STDELR-TF)
- 8" rear wheel 6" casters in front
- Back folds down for storage and transport and features deluxe back release (Figure C)
- Comes with aluminum rear wheel locks
- Convenient cup holder can easily attach to chair Sold separately (Model# STDS1040S) (Figure D)
- Weight Capacity: 300 lbs

PRODUCTS

Description

DFL19-BL	Blue Frame and Blue Plaid Upholstery, 1/cs
DFL19-BLK	Black Frame and Black/White/Red Plaid Upholstery, 1/cs
DFL19-RD	Burgundy Frame and Black Upholstery, 1/cs

SPECIFICATIONS

Description

Overall Width	22" (open) 9" (closed)
Overall Height	36"
Seat-To-Floor Height	18"
Back Height	18"
Seat Dimensions	19" (W) x 15.5" (D)
Overall Length w/ Riggings	33"
Weight Without Front Riggings	16.2 lbs
Weight With Front Riggings	18.9 lbs
Carton Shipping Weight	23.8 lbs.
Weight Limit	300 lbs.

ACCESSORIES

Description

STDS1040S	Universal Cup Holder, 1/ea
-----------	----------------------------