

Deluxe Fly-Weight Aluminum Transport Chair with Removable Casters



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



RETAIL BOX

- Weighs only 19 lbs
- Aluminum frame for lightweight strength and durability
- Easy folding and storage in closet or trunk (Figure A)
- Composite wheels are lightweight and maintenance free
- 8" wheels in front and rear easily remove without tools (Figure B - C)
- Deluxe back release (Figure D)
- Deluxe, all-aluminum rear wheel locks
- Comes with carry pocket on backrest (Figure E)
- Both 17" and 19" come with tool-free, swing-away footrests
- Comes with seat belt for added safety
- Flared arms allow the back to be folded virtually flat (Figure F)
- Soft, plush nylon upholstery
- Back folds down for easy storage and transport
- Convenient cup holder can easily attach to chair
- Sold separately (Model# STDS1040S)
- Weight Capacity: 300 lbs

PRODUCTS

PRODUCTS	Description
RTLFW19RW-RD	Red Frame, Black Upholstery, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Overall Width	22" (open) 9" (closed)
Overall Height	37.25"
Seat To-Floor Height	19"
Back Height	18"
Seat Dimensions	19" (W) x 16" (D)
Overall Length w/Riggings	33"
Weight	19 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES

ACCESSORIES	Description
STDS1040S	Universal Cup Holder, 1/ea