

22" Bariatric Aluminum Transport Chair With 12" Rear "Flat-Free" Wheels



FIGURE A



FIGURE B

- Weighs only 33 lbs
- Supports individuals up to 450 lbs
- Attachable push bar with padded foam grip (Figure A)
- 12" rear flat-free wheels
- Comes with seat belt
- Comes with swing-away footrests and anti-tippers
- Back folds down for convenient storage (Figure B)
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery
- Comes with hand wheel lock lever and companion activated wheel locks

PRODUCTS

Description

ATC22-R	Red, 1/cs
---------	-----------

SPECIFICATIONS

Description

Overall Width	30" (Open) 12" (Closed)
Overall Height	37.5"
Seat To-Floor Height	21"
Back Height	16"
Seat Dimensions	22" (W) x 18" (D)
Overall Length w/ Riggings	41"
Weight	33 lbs.
Carton Shipping Weight	39 lbs.
Weight Limit	450 lbs.
Warranty	Limited Lifetime