

Bariatric Sentra Extra-Heavy-Duty Wheelchair Dual Cross Brace

HCPCS: K0007



- Weight Capacity: 500 lbs
- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Reinforced side frame, caster journals and dual cross bars
- Dual-position front fork
- Dual-position front frame
- Composite Mag-style wheels with chrome hand rim are lightweight and maintenance free
- Hemi-adjustable caster is optional (Part #STDS1074)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Padded armrests provide added patient comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

SPECIFICATIONS	Description
Width Open with 20" Seat	29"
Width Open with 22" Seat	31"
Width Open with 24" Seat	32.5"
Width Closed	13"
Height	37"
Seat-To-Floor Height	17.5" - 19.5"
Back Height	18"
Seat Depth	18"
Length with Rigging	42"
Product Weight with 20" Seat	54 lbs.
Product Weight with 22" Seat	56 lbs.
Product Weight with 24" Seat	66.25 lbs.
Carton Shipping Weight with 20" Seat	60 lbs.
Carton Shipping Weight with 22" Seat	62 lbs.
Carton Shipping Weight with 24" Seat	72 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS801	Wheel Lock Extension, 1/ea
STDS802	Anti-Tippers With Wheels, 1/pr
STDS806	Anti-Fold Bar, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea

ACCESSORIES

Description

STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS2420	Seat Rail Extension Kit, 24" models only, 1/pr
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS862	Universal Forward Stabilizer, 1/pr
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDELR-AL	Articulating Legrest For Manual Wheelchairs, 1/cs

PARTS

Description

STDS801	Wheel Lock Extension, 1/ea
STDS802	Anti-Tippers With Wheels, 1/pr
STDS806	Anti-Fold Bar, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS836	Chart Carry Pocket, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS2420	Seat Rail Extension Kit, 24" models only, 1/pr
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS862	Universal Forward Stabilizer, 1/pr
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS2329BLK	1/bx
STDS2331BLK	1/bx
STDS2333BLK	1/bx
STDS871	1/ea
STDS1000HD	1/ea
STDS1052	1/ea
STDS1012HD	1/ea
STDS005008	1/ea
STDS4020	1/ea
STDS2A5222	1/ea
STDS001061	1/ea
STDS4011	1/ea
STDS4010R	1/ea
SIDS4010L	
STDS1002-1	1/ea
STDS1011	1/ea
STDS1075	1/ea
STDS1074	1/ea