



Chrome Sport Wheelchair

Dual axle

HCPCS: K0001 / K0002 / K0007



INSET A

- Carbon steel frame with chrome coating provides an attractive, chip-proof, easy-to-maintain finish (Figure A)
- Embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Urethane tires mounted on composite wheels provide durability and low maintenance
- 8" front casters, adjustable in three positions
- Padded armrests provide additional comfort
- Carry pocket on backrest provides additional convenience
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Plastic foot plates
- Comes with push-to-lock wheel locks

PRODUCTS

Description

CS16DDA-SF	ChromeSport W/C 16Det. Desk SF, 1/cs
CS16DDA-ELR	ChromeSport W/C 16Det. Desk ELR, 1/cs
CS18DDA-SF	Chrome Sport W/C 18 DET.DESK, 1/cs
CS18DDA-ELR	Chrome Sport W/C 18 DET. DESK ELR, 1/cs
CS20DDA-SF	ChromeSport W/C 20in Det. Desk SF, 1/cs
CS20DDA-ELR	ChromeSport W/C 20Det. Desk ELR, 1/cs
CS20DFA-SF	ChromeSport W/C 20Det.FullArm SF, 1/cs
CS20DFA-ELR	ChromeSport W/C 20Det.FullArm ELR, 1/cs

SPECIFICATIONS

Description

Width Open with 16" Seat	24"
Width Open with 18" Seat	26"
Width Open with 20" Seat	28"
Width Closed	12.5"
Height	36"
Seat-To-Floor Height	17.5" - 19.5"
Back Height	16"
Seat Depth	16"
Length with Rigging	42"
Product Weight with 16" Seat	39 lbs.
Product Weight with 18" Seat	40 lbs.
Product Weight with 20" Seat	41 lbs.
Carton Shipping Weight with 16" Seat	45 lbs.
Carton Shipping Weight with 18" Seat	46 lbs.
Carton Shipping Weight with 20" Seat	47 lbs.
Weight Limit with 16" Seat	250 lbs.
Weight Limit with 18" Seat	300 lbs.

SPECIFICATIONS	Description
Weight Limit with 20" Seat	350 lbs.
Warranty	Limited 3 Year

ACCESSORIES	Description
STDS806	Anti-Fold Bar, 1/ea
STDS833N	Anti-Tippers With Wheels Black, 1/pr
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea

PARTS	Description
STDS806	Anti-Fold Bar, 1/ea
STDS833N	Anti-Tippers With Wheels Black, 1/pr
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS836	Chart Carry Pocket, 1/ea
STDS831	Heel Loops, 1/bx
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea