



## Deluxe Sentra Full Reclining Wheelchair Dual Axle Wheelchair

HCPCS: K0001 / E1226

- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Durable, heavy-gauge naugahyde upholstery
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Composite Mag-style wheels with chrome hand rims are lightweight and maintenance free
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels set back on frame prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers standard
- Comes standard with swing-away elevating legrests
- Aluminum foot plates are attractive, crack-proof and lightweight
- Front caster forks are adjustable in two positions
- Carry pocket standard
- Headrest extension with cushioned head immobilizer standard
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD22RBDDA	Sentra 22Recl.DetDskArm, 1/cs

SPECIFICATIONS	Description
Width Open with 22" Seat	30"
Width Closed (with 14", 16", 18", 20", 22")	12.5"
Height with 22" Seat	51"
Seat-To-Floor Height	21"
Back Height	23"   33"
Length with Riggings	51"
Product Weight with 22" Seat	62 lbs.
Carton Shipping Weight with 22" Seat	72 lbs.
Weight Limit with 22" Seat	450 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS801	Wheel Lock Extension, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS802	Anti-Tippers With Wheels, 1/pr
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS840	Flat-free Wheel With Handrims, 1/ea
STDS831	Heel Loops, 1/bx
STDS1001	Standard 24" Wheel, 1/ea

## ACCESSORIES

## Description

WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS1022	Oblique Projection Handrim, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDSC1818V	Seat Rail Extension Kit, 1/pr
STDS1002	Standard 8" Caster, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS6008-1	Oxygen Cylinder Carry Bag for Wheelchair, 1/ea
STDS4022R	Wheel Locks, Right, Pull to Lock, Side Mount, Detachable Arm, 1/ea
STDS4022L	Wheel Locks, Left, Pull to Lock, Side Mount, Detachable Arm, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS4035R	Wheel Locks, Right, Push to Lock, Side Mount, Detachable Arm, 1/ea
STDS4035L	Wheel Locks, Left, Push to Lock, Side Mount, Detachable Arm, 1/ea
STDS1074	6" Caster Wheel, 1/ea

## PARTS

## Description

STDS801	Wheel Lock Extension, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS802	Anti-Tippers With Wheels, 1/pr
STDS818	Anti-Tippers without Wheels, Chrome, 1/pr
STDS835	1/ea
STDS840	Flat-free Wheel With Handrims, 1/ea
STDS831	Heel Loops, 1/bx
STDS1001	Standard 24" Wheel, 1/ea
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS1022	Oblique Projection Handrim, 1/ea
STDS823	Overhead Anti-Theft Device, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS811	Seat Rail Extensions, 1/pr
STDS1002	Standard 8" Caster, 1/ea
STDS860	Universal Forward Stabilizer, 1/pr
STDS6008-1	Oxygen Cylinder Carry Bag for Wheelchair, 1/ea
STDS4022R	Wheel Locks, Right, Pull to Lock, Side Mount, Detachable Arm, 1/ea
STDS4022L	Wheel Locks, Left, Pull to Lock, Side Mount, Detachable Arm, 1/ea
STDS4035R	Wheel Locks, Right, Push to Lock, Side Mount, Detachable Arm, 1/ea
STDS4035L	Wheel Locks, Left, Push to Lock, Side Mount, Detachable Arm, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea