

Specifically Designed for Mid-Wheel Power Chairs

Compact, open platform lift is designed to easily transport 10" or 14" mid-wheel drive (MWD) power chairs with solid suspension casters. Handles capacities up to 350 lbs.

Features & Benefits

- Padded hold-down arm conveniently secures power wheelchair
- No chair modifications required
- Incredibly small platform, just 16.25" x 27"
- Fully powered
- Class II option available
- Swing-Away option available



Specifications	AL580	AL580XL Extra Large
Accommodates	Mid-wheel power chairs with 10" drive wheels	Mid-wheel power chairs with 14" drive wheels
Platform	16.5" x 27"	17.5" x 27"
Lifting Capacity	350 lbs. (160 kg)	350 lbs. (160 kg)
Installed Weight	63 lbs. (31 kg)	63 lbs. (31 kg)
Securement	Padded hold-down arm & wheel cradle	Padded hold-down arm & wheel cradle
Hitch	Class II or greater	Class II or greater
Included Accessories	AL123 Class III hitch adapter / AL107 hitch post	AL123 Class III hitch adapter / AL107 hitch post
Options	Battery pack, Swing-Away, Power chair covers	Battery pack, Swing-Away, Power chair covers
Warranty	3-year transferable	3-year transferable